MEDICAL DEVICE REPROCESSING ASSOCIATION OF ONTARIO (MDRAO) PRACTICE EXAM (SAMPLE)

STUDY GUIDE



SAMPLE STUDY GUIDE WITH LIMITED QUESTIONS

VISIT US ONLINE FOR THE FULL VERSION STUDY GUIDE

https://mdrao.examzify.com

Copyright © 2026 by Examzify - A Kaluba Technologies Inc. product.

ALL RIGHTS RESERVED.

No part of this book may be reproduced or transferred in any form or by any means, graphic, electronic, or mechanical, including photocopying, recording, web distribution, taping, or by any information storage retrieval system, without the written permission of the author.

Notice: Examzify makes every reasonable effort to obtain accurate, complete, and timely information about this product from reliable sources.

SAMPLE

Table of Contents

Copyright 1

Table of Contents 2

Introduction 3

How to Use This Guide 4

Questions 5

Answers 8

Explanations 10

Next Steps 15

SAMPLE

Introduction

Preparing for a certification exam can feel overwhelming, but with the right tools, it becomes an opportunity to build confidence, sharpen your skills, and move one step closer to your goals. At Examzify, we believe that effective exam preparation isn't just about memorization, it's about understanding the material, identifying knowledge gaps, and building the test-taking strategies that lead to success.

This guide was designed to help you do exactly that.

Whether you're preparing for a licensing exam, professional certification, or entry-level qualification, this book offers structured practice to reinforce key concepts. You'll find a wide range of multiple-choice questions, each followed by clear explanations to help you understand not just the right answer, but why it's correct.

The content in this guide is based on real-world exam objectives and aligned with the types of questions and topics commonly found on official tests. It's ideal for learners who want to:

- Practice answering questions under realistic conditions,
- Improve accuracy and speed,
- Review explanations to strengthen weak areas, and
- Approach the exam with greater confidence.

We recommend using this book not as a stand-alone study tool, but alongside other resources like flashcards, textbooks, or hands-on training. For best results, we recommend working through each question, reflecting on the explanation provided, and revisiting the topics that challenge you most.

Remember: successful test preparation isn't about getting every question right the first time, it's about learning from your mistakes and improving over time. Stay focused, trust the process, and know that every page you turn brings you closer to success.

Let's begin.

SAMPLE

How to Use This Guide

This guide is designed to help you study more effectively and approach your exam with confidence. Whether you're reviewing for the first time or doing a final refresh, here's how to get the most out of your Examzify study guide:

1. Start with a Diagnostic Review

Skim through the questions to get a sense of what you know and what you need to focus on. Your goal is to identify knowledge gaps early.

2. Study in Short, Focused Sessions

Break your study time into manageable blocks (e.g. 30 – 45 minutes). Review a handful of questions, reflect on the explanations.

3. Learn from the Explanations

After answering a question, always read the explanation, even if you got it right. It reinforces key points, corrects misunderstandings, and teaches subtle distinctions between similar answers.

4. Track Your Progress

Use bookmarks or notes (if reading digitally) to mark difficult questions. Revisit these regularly and track improvements over time.

5. Simulate the Real Exam

Once you're comfortable, try taking a full set of questions without pausing. Set a timer and simulate test-day conditions to build confidence and time management skills.

6. Repeat and Review

Don't just study once, repetition builds retention. Re-attempt questions after a few days and revisit explanations to reinforce learning. Pair this guide with other Examzify tools like flashcards, and digital practice tests to strengthen your preparation across formats.

There's no single right way to study, but consistent, thoughtful effort always wins. Use this guide flexibly, adapt the tips above to fit your pace and learning style. You've got this!

SAMPLE

Questions

SAMPLE

1. Which practice is essential for effective hand hygiene?

A. Using gloves only

B. Washing hands with soap and water

C. Applying lotion afterward

D. Using disinfectants exclusively

2. What distinguishes soft water from hard water?

A. Soft water is mineralized

B. Hard water is pure

C. Soft water lacks minerals

D. Hard water contains only magnesium

3. Can dry scrubbers be used to clean areas in sterile storage?

A. Yes, they are effective

B. No, only wet scrubbers may be used

C. Only for non-critical areas

D. Only if they are HEPA-filtered

4. The steam trap in a sterilizer is primarily responsible for what?

A. Regulating steam intake

B. Removing air at the beginning

C. Monitoring temperature cycles

D. Conserving energy

5. What is the purpose of a laryngoscope blade?

A. To examine the nose

B. To examine the ears

C. To examine the throat

D. To perform surgery

6. What is the humidity range required for storage temperatures in Tier 1?

A. 20-40%

B. 30-50%

C. 30-60%

D. 40-70%

SAMPLE

7. What does the Bowie-Dick Type Test assess?

A. The integrity of sterilization packaging

B. The presence of air in the sterilization package

C. The effectiveness of cleaning agents

D. The temperature accuracy of the sterilizer

8. How should reusable brushes be handled after use?

A. Discarded immediately

B. Cleaned and disinfected after each use

C. Rinsed with water only

D. Stored dry for future use

9. How does steam sterilization differ from dry heat sterilization in terms of required temperatures?

A. Requires lower temperatures for shorter times

B. Requires the same temperatures for longer times

C. Requires higher temperatures for shorter times

D. Requires higher temperatures for longer times

10. What are the two types of sterile packaging wrapping?

A. Envelope and square wrapping

B. Flat and tubular wrapping

C. Plastic and cloth wrapping

D. Open and closed wrapping

SAMPLE

Answers

SAMPLE

1. B
2. C
3. B
4. A
5. C
6. C
7. B
8. B
9. A
10. A

SAMPLE

Explanations

SAMPLE

1. Which practice is essential for effective hand hygiene?

A. Using gloves only

B. Washing hands with soap and water

C. Applying lotion afterward

D. Using disinfectants exclusively

Washing hands with soap and water is essential for effective hand hygiene because it is the method most effective in removing dirt, bacteria, viruses, and other contaminants from the skin. This practice is supported by health organizations worldwide, as it significantly reduces the likelihood of transmitting infections. Soap works by breaking down oils and dirt, allowing them to be rinsed away with water, which provides a mechanical action that further aids in removing pathogens. While using gloves can add a layer of protection and applying lotion afterward can help maintain skin integrity, neither are substitutes for handwashing. Disinfectants can also be effective under certain circumstances, but their proper use may not achieve the same thoroughness in removing organic materials as washing with soap and water does. Therefore, regular handwashing remains a critical component of infection control practices.

2. What distinguishes soft water from hard water?

A. Soft water is mineralized

B. Hard water is pure

C. Soft water lacks minerals

D. Hard water contains only magnesium

Soft water is characterized by its low mineral content, specifically the absence or low levels of calcium and magnesium ions, which are the minerals that contribute to water hardness. The distinction between soft and hard water primarily revolves around these mineral concentrations. Soft water results from processes that remove these minerals, such as ion exchange or filtration, making it gentler on skin and appliances, and it tends to lather more easily with soap. This is why the definition of soft water focuses on its lack of minerals, as opposed to being mineralized or being "pure" in a different sense. In contrast, hard water is defined by its higher concentration of minerals, particularly calcium and magnesium, which can cause scaling in pipes and reduction in soap effectiveness. The presence of only magnesium in hard water is not accurate, as hard water can contain various other minerals as well, leading to variability in hardness levels. Thus, the correct answer consistently highlights the main characteristic of soft water being its lack of significant mineral content.

3. Can dry scrubbers be used to clean areas in sterile storage?

A. Yes, they are effective

B. No, only wet scrubbers may be used

C. Only for non-critical areas

D. Only if they are HEPA-filtered

Dry scrubbers are not appropriate for cleaning areas in sterile storage because they can potentially introduce contaminants, such as dust or particulate matter, into the sterile environment. Sterile storage areas require a controlled environment where maintaining sterility is paramount. Using dry scrubbers could disturb settled particles and could inadvertently compromise the sterility of the items stored. In contrast, wet scrubbers are effective for cleaning as they typically use a liquid solution, which can help to reduce airborne particles and promote a more thorough cleaning process that maintains the sterile environment. Wet cleaning methods tend to minimize the risk of contamination by ensuring surfaces are cleaner and less likely to harbor particles that could compromise sterility. Thus, understanding the appropriate methods of cleaning sterile storage areas is critical for infection control and overall patient safety, highlighting why the use of dry scrubbers is not recommended.

SAMPLE

4. The steam trap in a sterilizer is primarily responsible for what?

A. Regulating steam intake

B. Removing air at the beginning

C. Monitoring temperature cycles

D. Conserving energy

The steam trap in a sterilizer primarily functions to remove condensate and non-condensable gases, such as air, from the sterilization chamber. This role is crucial for maintaining optimal conditions within the sterilizer, as the presence of air can interfere with the penetration of steam and ultimately impact the effectiveness of the sterilization process. In this context, one of the most critical functions of the steam trap is to ensure that air is removed at the beginning of the cycle, allowing for proper steam penetration and heat transfer to the items being sterilized. By efficiently removing air, the steam trap helps in achieving the desired temperature and pressure necessary for effective sterilization. Thus, while the regulation of steam intake and energy conservation are important, they are secondary to the immediate need for air removal at the start of the sterilization process, which is essential for ensuring that the sterilization parameters are met. This highlights the steam trap's primary responsibility in facilitating an effective and reliable sterilization process, which includes the timely removal of air from the chamber to allow for maximum steam efficiency.

5. What is the purpose of a laryngoscope blade?

A. To examine the nose

B. To examine the ears

C. To examine the throat

D. To perform surgery

A laryngoscope blade is specifically designed for examining the throat, particularly to visualize the larynx and vocal cords. It serves as an essential tool in procedures such as intubation and in assessing airway patency, making it indispensable in both emergency and routine medical situations involving the upper airway. The structure of the laryngoscope allows healthcare providers to lift the tongue and expand the view of the throat, providing a clear line of sight to assist in various medical interventions. This specificity in function is what distinguishes it from tools intended for examining the nose or ears, or for performing surgeries, which would employ different types of instruments tailored for those activities.

6. What is the humidity range required for storage temperatures in Tier 1?

A. 20-40%

B. 30-50%

C. 30-60%

D. 40-70%

The correct answer for the humidity range required for storage temperatures in Tier 1 is 30-60%. This range is specified to ensure that stored medical devices maintain their integrity, functionality, and sterility over time. Proper humidity control in storage environments is crucial because excessive moisture can lead to microbial growth, while too little humidity can cause some materials to degrade or become brittle. In the context of medical device reprocessing, maintaining the appropriate humidity levels helps prevent damage to packaging and assists in preserving the sterility of the devices until they are needed for use, ensuring patient safety. The 30-60% range has been established based on best practices and guidelines within the industry to balance these concerns effectively. Keeping humidity within this range supports the overall aim of the MDRAO to promote safe and effective storage practices for medical devices.

SAMPLE

7. What does the Bowie-Dick Type Test assess?

A. The integrity of sterilization packaging

B. The presence of air in the sterilization package

C. The effectiveness of cleaning agents

D. The temperature accuracy of the sterilizer

The Bowie-Dick Type Test is specifically designed to assess the presence of air in the sterilization chamber during steam sterilization processes. Steam sterilization relies on the effective penetration of steam to ensure that all surfaces of the items being sterilized come into contact with the steam at the appropriate temperature and pressure. If air is present in the sterilization chamber, it can create barriers that prevent proper steam penetration, compromising the sterilization cycle and resulting in ineffective sterilization. The test utilizes a special test pack that contains a porous medium and a chemical indicator. If air is present, it will cause the test pack to fail, indicated by the reactive indicator not changing color. Therefore, detecting air in the sterilization chamber is critical for confirming that the sterilization process will be effective. While other options touch upon relevant concepts in sterilization, they do not align with the specific purpose of the Bowie-Dick Type Test. For example, assessing the integrity of sterilization packaging involves checking for physical damage or leaks, but not the presence of air in the chamber itself. The effectiveness of cleaning agents is related to the cleaning process before sterilization, and temperature accuracy ensures that the sterilizer operates within specified parameters, but these are not the focus of the Bowie-Dick Type

8. How should reusable brushes be handled after use?

A. Discarded immediately

B. Cleaned and disinfected after each use

C. Rinsed with water only

D. Stored dry for future use

Reusable brushes can carry microbes from one instrument to the next, so after each use they need to be cleaned to remove organic material and then disinfected to kill remaining microorganisms. This two-step approach reduces bioburden and prevents cross-contamination, protecting patients and staff. The proper sequence is to rinse to remove loose debris, wash with an appropriate detergent, rinse thoroughly, then apply a validated disinfectant with the correct contact time according to policy and manufacturer instructions, and finally dry completely before storage. Rinsing with water alone or storing a used brush without disinfection would leave contaminants and promote microbial growth, and discarding a brush after every use is wasteful and unnecessary when cleaning and disinfection are effective.

SAMPLE

9. How does steam sterilization differ from dry heat sterilization in terms of required temperatures?

A. Requires lower temperatures for shorter times

B. Requires the same temperatures for longer times

C. Requires higher temperatures for shorter times

D. Requires higher temperatures for longer times

Steam sterilization differs from dry heat sterilization in that it requires lower temperatures while achieving effective sterilization in shorter times. This is primarily due to the presence of moisture in steam, which significantly enhances its heat transfer properties. The steam can penetrate the materials being sterilized more effectively than hot air, allowing it to achieve the same sterilization outcomes at lower temperatures. In contrast, dry heat sterilization typically operates at higher temperatures for longer durations to achieve the same sterilization effect. The lack of moisture means that more energy is required to kill microorganisms, resulting in the necessity for elevated temperatures and extended exposure times. Understanding this difference is crucial for selecting the appropriate sterilization method for various medical devices while ensuring compliance with safety standards and effective infection control practices.

10. What are the two types of sterile packaging wrapping?

A. Envelope and square wrapping

B. Flat and tubular wrapping

C. Plastic and cloth wrapping

D. Open and closed wrapping

The two types of sterile packaging wrapping commonly recognized in the field of medical device reprocessing are flat wrapping and tubular wrapping. Flat wrapping refers to the technique where items are wrapped flatly in a specific material, allowing for easy placement in sterilization equipment, providing reliable sealing and sterilization. Tubular wrapping typically involves wrapping items that are long or cylindrical, accommodating for their shape and ensuring that they remain sterile until they are needed. While envelope and square wrapping may be understood in a similar context, they do not accurately represent the commonly accepted terminology used in reprocessing practices. It's important to focus on the recognized categories—flat and tubular—as they reflect established methods in sterilization processes used to maintain the sterility of medical devices. This terminology enhances effective communication in the field and supports proper training in sterile processing techniques.

SAMPLE

Next Steps

Congratulations on reaching the final section of this guide. You've taken a meaningful step toward passing your certification exam and advancing your career.

As you continue preparing, remember that consistent practice, review, and self-reflection are key to success. Make time to revisit difficult topics, simulate exam conditions, and track your progress along the way.

If you need help, have suggestions, or want to share feedback, we'd love to hear from you. Reach out to our team at hello@examzify.com.

Or visit your dedicated course page for more study tools and resources:

https://mdrao.examzify.com

We wish you the very best on your exam journey. You've got this!

SAMPLE