# REGISTERED SLEEP TECHNOLOGIST PRACTICE EXAM (SAMPLE)

## STUDY GUIDE



## SAMPLE STUDY GUIDE WITH LIMITED QUESTIONS

**VISIT US ONLINE FOR THE FULL VERSION STUDY GUIDE**

https://registeredsleeptechnologist.examzify.com

**Copyright © 2026 by Examzify - A Kaluba Technologies Inc. product.**

**ALL RIGHTS RESERVED.**

No part of this book may be reproduced or transferred in any form or by any means, graphic, electronic, or mechanical, including photocopying, recording, web distribution, taping, or by any information storage retrieval system, without the written permission of the author.

Notice: Examzify makes every reasonable effort to obtain accurate, complete, and timely information about this product from reliable sources.

SAMPLE

# Table of Contents


Copyright ..... 1

Table of Contents ..... 2

Introduction ..... 3

How to Use This Guide ..... 4

Questions ..... 5

Answers ..... 8

Explanations ..... 10

Next Steps ..... 16


SAMPLE

# Introduction

Preparing for a certification exam can feel overwhelming, but with the right tools, it becomes an opportunity to build confidence, sharpen your skills, and move one step closer to your goals. At Examzify, we believe that effective exam preparation isn't just about memorization, it's about understanding the material, identifying knowledge gaps, and building the test-taking strategies that lead to success.

This guide was designed to help you do exactly that.

Whether you're preparing for a licensing exam, professional certification, or entry-level qualification, this book offers structured practice to reinforce key concepts. You'll find a wide range of multiple-choice questions, each followed by clear explanations to help you understand not just the right answer, but why it's correct.

The content in this guide is based on real-world exam objectives and aligned with the types of questions and topics commonly found on official tests. It's ideal for learners who want to:

- Practice answering questions under realistic conditions,
- Improve accuracy and speed,
- Review explanations to strengthen weak areas, and
- Approach the exam with greater confidence.

We recommend using this book not as a stand-alone study tool, but alongside other resources like flashcards, textbooks, or hands-on training. For best results, we recommend working through each question, reflecting on the explanation provided, and revisiting the topics that challenge you most.

Remember: successful test preparation isn't about getting every question right the first time, it's about learning from your mistakes and improving over time. Stay focused, trust the process, and know that every page you turn brings you closer to success.

Let's begin.

SAMPLE

# How to Use This Guide

This guide is designed to help you study more effectively and approach your exam with confidence. Whether you're reviewing for the first time or doing a final refresh, here's how to get the most out of your Examzify study guide:

### 1. Start with a Diagnostic Review

Skim through the questions to get a sense of what you know and what you need to focus on. Your goal is to identify knowledge gaps early.

### 2. Study in Short, Focused Sessions

Break your study time into manageable blocks (e.g. 30 – 45 minutes). Review a handful of questions, reflect on the explanations.

### 3. Learn from the Explanations

After answering a question, always read the explanation, even if you got it right. It reinforces key points, corrects misunderstandings, and teaches subtle distinctions between similar answers.

### 4. Track Your Progress

Use bookmarks or notes (if reading digitally) to mark difficult questions. Revisit these regularly and track improvements over time.

### 5. Simulate the Real Exam

Once you're comfortable, try taking a full set of questions without pausing. Set a timer and simulate test-day conditions to build confidence and time management skills.

### 6. Repeat and Review

Don't just study once, repetition builds retention. Re-attempt questions after a few days and revisit explanations to reinforce learning. Pair this guide with other Examzify tools like flashcards, and digital practice tests to strengthen your preparation across formats.

There's no single right way to study, but consistent, thoughtful effort always wins. Use this guide flexibly, adapt the tips above to fit your pace and learning style. You've got this!

SAMPLE

# Questions

SAMPLE

**1. What is the frequency range of alpha waves?**

A. 1-3 Hz

B. 4-7 Hz

C. 8-13 Hz

D. 14-30 Hz

**2. Which of the following is an example of a direct current channel?**

A. EEG

B. SpO2

C. EKG

D. EMG

**3. What does the recording of eye movements in sleep studies measure?**

A. The electrical potential difference between the cornea and the retina

B. The movement of eyelids

C. The frequency of eye blinks

D. Changes in pupil size

**4. What is the primary difference between PLMD and RLS?**

A. PLMD occurs during sleep, while RLS can occur during wakefulness

B. RLS is sleep-related and PLMD is not

C. RLS is always accompanied by pain

D. PLMD can be treated with medications but RLS cannot

**5. How do patients with a hypoxic drive typically react to CPAP?**

A. They respond well to CPAP

B. They do not respond well to CPAP

C. They require lower CPAP settings

D. They are unaffected by CPAP

**6. What physiological state is typically associated with alpha waves in an EEG?**

A. Deep sleep

B. Active concentration

C. Relaxed wakefulness

D. Rapid eye movement sleep

SAMPLE

**7. What is the term used when electrical signals from external sources interfere with signals from the patient?**

A. Stray capacitance

B. Noise artifact

C. Signal distortion

D. Electromagnetic interference

**8. Sleep hygiene practices primarily aim to promote which of the following?**

A. Short-term sleep deprivation

B. Consistent sleep routines

C. Medication dependence

D. Extended wakefulness

**9. Which of the following is NOT considered a predisposing factor for Obstructive Sleep Apnea (OSA)?**

A. Obesity

B. Epileptiform activity

C. Aging

D. Alcohol use

**10. What filter settings are most appropriate for EMG channels?**

A. LFF=5 Hz HFF=50 Hz

B. LFF=10 Hz HFF=100 Hz

C. LFF=20 Hz HFF=1000 Hz

D. LFF=1 Hz HFF=30 Hz

SAMPLE

# Answers

SAMPLE

1. C
2. B
3. A
4. A
5. B
6. C
7. A
8. B
9. B
10. B

SAMPLE

# Explanations

SAMPLE

## 1. What is the frequency range of alpha waves?

A. 1-3 Hz

B. 4-7 Hz

C. 8-13 Hz

D. 14-30 Hz

*Alpha waves are a specific type of brainwave pattern that are typically associated with a relaxed, yet alert state of mind. They are most prominent when a person is awake but relaxed, often occurring during moments of calmness like meditation or light rest. The frequency range of alpha waves is classified between 8 and 13 Hz. This range is distinctive because it lies at a transitional frequency between the slower theta waves and the faster beta waves. Understanding this frequency is crucial for sleep technologists, as alpha waves are indicative of certain cognitive states and can have implications for sleep staging and the identification of sleep disorders. Recognizing alpha waves during EEG readings can help in interpreting the mental state of the subject being studied. In this context, the other frequency ranges provided do not align with the characteristics of alpha waves. The ranges listed for 1-3 Hz correspond to delta waves, which are linked to deep sleep, while 4-7 Hz pertains to theta waves, typically associated with light sleep and relaxation. On the other end, 14-30 Hz relates to beta waves, commonly observed during active, busy, or anxious states of mind. Therefore, knowing that alpha waves fall specifically within the 8-13 Hz range allows professionals to accurately assess brain*

## 2. Which of the following is an example of a direct current channel?

A. EEG

B. SpO2

C. EKG

D. EMG

*The correct choice identifies SpO2 as an example of a direct current channel because it measures the oxygen saturation level in the blood using a photoplethysmography technique, which relies on continuous waveforms produced by the absorption of infrared light. This type of measurement is consistent and does not involve alternating currents, hence qualifying it as a direct current. In contrast, the other options illustrate different modalities. EEG (electroencephalography) captures electrical activity in the brain and typically utilizes alternating current to assess brain wave patterns. EKG (electrocardiography) measures the electrical activity of the heart and often involves alternating current signals as heart rhythms fluctuate. EMG (electromyography) measures the electrical activity produced by skeletal muscles and also works with signals that can vary over time, characteristic of alternating current systems. Understanding these distinctions clarifies why SpO2 is the only direct current example among the provided choices.*

SAMPLE

## 3. What does the recording of eye movements in sleep studies measure?

A. The electrical potential difference between the cornea and the retina

B. The movement of eyelids

C. The frequency of eye blinks

D. Changes in pupil size

*The recording of eye movements in sleep studies primarily measures the electrical potential difference between the cornea and the retina. This phenomenon is referred to as electrooculography (EOG). In this process, electrodes placed around the eyes detect the electrical signals generated by the different charges within the eye. Since the cornea has a positive charge relative to the retina, the movement of the eyes creates a change in electrical potential that can be captured for analysis. Understanding eye movement is crucial during sleep studies, particularly in the identification of different sleep stages. Rapid Eye Movement (REM) sleep, characterized by its unique brain wave patterns and distinct eye movements, can be distinguished using this method. The ability to assess these eye movements contributes valuable information to the interpretation of various sleep phenomena, including sleep disorders. Other options, such as the movement of eyelids, frequency of eye blinks, or changes in pupil size, do not directly correlate with the measurements generally taken during sleep studies. While these actions may occur during sleep, they do not provide the same level of insight into sleep stages as the electrical potential difference noted in electrooculography. Thus, measuring the electrical potential difference is the most meaningful and relevant approach in the context of sleep studies.*

## 4. What is the primary difference between PLMD and RLS?

A. PLMD occurs during sleep, while RLS can occur during wakefulness

B. RLS is sleep-related and PLMD is not

C. RLS is always accompanied by pain

D. PLMD can be treated with medications but RLS cannot

*The primary difference between Periodic Limb Movement Disorder (PLMD) and Restless Legs Syndrome (RLS) lies in when these conditions manifest. PLMD is characterized by repetitive movements of the limbs during sleep, and these movements often occur in clusters, leading to disrupted sleep. In contrast, RLS involves an overwhelming urge to move the legs, which typically arises in periods of rest or inactivity, and it can be experienced while awake. This distinction highlights that while PLMD specifically occurs during sleep, RLS can be triggered at any time when an individual is awake, especially during times of sedentary activity. This feature of RLS makes it uniquely tied to the sensations of discomfort that prompt the need to move, providing a clear differentiation in symptom presentation between the two disorders. In summary, the recognition that PLMD is confined to sleep while RLS can manifest during wakefulness is crucial for proper diagnosis and treatment strategies.*

SAMPLE

## 5. How do patients with a hypoxic drive typically react to CPAP?

A. They respond well to CPAP

B. They do not respond well to CPAP

C. They require lower CPAP settings

D. They are unaffected by CPAP

*Patients with a hypoxic drive typically do not respond well to CPAP. The hypoxic drive refers to the physiological response where individuals rely on low oxygen levels (rather than high carbon dioxide levels) to stimulate breathing. In conditions such as chronic obstructive pulmonary disease (COPD), these patients may have a depressed central respiratory drive, meaning that increasing airway pressure with CPAP could reduce their urge to breathe because it is counteracting their natural response to low oxygen levels. When CPAP is applied, it can increase the overall oxygen level, which may lead to a further decrease in the respiratory drive. This can cause an inadequate response to the treatment, making it more challenging for these patients to achieve optimal ventilation and oxygenation. Therefore, understanding the dynamics of the hypoxic drive is crucial when managing such patients with CPAP, as it can have adverse effects rather than providing the intended support.*

## 6. What physiological state is typically associated with alpha waves in an EEG?

A. Deep sleep

B. Active concentration

C. Relaxed wakefulness

D. Rapid eye movement sleep

*Alpha waves in an EEG reflect a specific brain state that occurs during relaxed wakefulness. When an individual is awake but in a calm, restful state, such as when they are meditating or calmly resting with their eyes closed, alpha waves are prominent. This EEG pattern indicates a decrease in mental activity compared to fully alert states, where beta waves dominate. In contrast, deep sleep is characterized by delta waves, which signify a much lower level of consciousness and brain activity. Active concentration generally correlates with beta waves, which are associated with alertness and intense mental activity. Rapid eye movement (REM) sleep, on the other hand, is associated with beta waves and sawtooth waves and occurs after cycles of deep sleep. Thus, the physiological state associated with alpha waves is indeed relaxed wakefulness, as this state represents a balancced transition towards deeper sleep or reductions in cognitive function without the individual being fully asleep.*

SAMPLE

## 7. What is the term used when electrical signals from external sources interfere with signals from the patient?

A. Stray capacitance

B. Noise artifact

C. Signal distortion

D. Electromagnetic interference

*The term that describes when electrical signals from external sources interfere with signals from the patient is noise artifact. This phenomenon occurs when unwanted signals are present in the data collected during a sleep study, which can obscure the true physiological signals that technologists are trying to measure. Noise artifacts can arise from a variety of sources, including electronic equipment, power lines, and even movement from the patient. Understanding noise artifacts is crucial for sleep technologists, as they need to identify and minimize these interferences to ensure accurate data collection and interpretation. Recognizing these artifacts helps in diagnosing sleep disorders more effectively. While stray capacitance refers to unintended capacitive coupling, and electromagnetic interference is more broadly related to disruptions caused by electromagnetic fields, noise artifacts specifically highlight the effect on the quality of the signals being analyzed in a clinical context. Signal distortion, on the other hand, is a more general term that can refer to any alteration of a signal, but it does not capture the specific nature of external interference affecting the collected data.*

## 8. Sleep hygiene practices primarily aim to promote which of the following?

A. Short-term sleep deprivation

B. Consistent sleep routines

C. Medication dependence

D. Extended wakefulness

*Sleep hygiene practices primarily aim to promote consistent sleep routines. This approach is designed to facilitate better sleep quality and enhance overall health by encouraging individuals to follow predictable sleep patterns. Establishing regular sleep schedules helps regulate the body's internal clock, known as the circadian rhythm, which plays a vital role in the sleep-wake cycle. By promoting consistent sleep routines, individuals are more likely to fall asleep and wake up at similar times each day, which can lead to improved sleep efficiency and duration. Additionally, adhering to healthy sleep habits, such as creating a relaxing bedtime environment and avoiding stimulants before bed, further supports the effectiveness of these routines. The other options focus on aspects that are contrary to the goals of sleep hygiene. Short-term sleep deprivation can lead to various negative health consequences and is not a goal of sleep hygiene. Medication dependence, while sometimes a necessary intervention, is not something that sleep hygiene aims to encourage; rather, it seeks alternative methods to improve sleep. Finally, extended wakefulness is detrimental to health and well-being and stands in opposition to the recommendations made through proper sleep hygiene practices.*

SAMPLE

## 9. Which of the following is NOT considered a predisposing factor for Obstructive Sleep Apnea (OSA)?

A. Obesity

B. Epileptiform activity

C. Aging

D. Alcohol use

*Epileptiform activity is not considered a predisposing factor for Obstructive Sleep Apnea (OSA). OSA is primarily characterized by repeated episodes of upper airway obstruction during sleep, and its risk factors predominantly include conditions and behaviors that influence airway patency and sleep quality. Obesity is a significant risk factor due to the accumulation of adipose tissue around the neck, which can lead to a narrowing of the airway. Aging is another well-known factor, as physiological changes that occur with age, including muscle tone reduction, can contribute to increased resistance in the airway during sleep. Alcohol use can also exacerbate OSA by relaxing the throat muscles, leading to a higher likelihood of airway obstruction. In contrast, while epileptiform activity may lead to different sleep disturbances or conditions, it is not directly associated with the anatomical and physiological factors that contribute to the prevalence of OSA. This distinction highlights the importance of understanding which conditions and lifestyle factors are directly linked to the risk of obstructive sleep apnea.*

## 10. What filter settings are most appropriate for EMG channels?

A. LFF=5 Hz HFF=50 Hz

B. LFF=10 Hz HFF=100 Hz

C. LFF=20 Hz HFF=1000 Hz

D. LFF=1 Hz HFF=30 Hz

*The most appropriate filter settings for EMG (electromyography) channels are 10 Hz for the low frequency filter (LFF) and 100 Hz for the high frequency filter (HFF). These settings are chosen to effectively capture the relevant frequency components of muscle activity while minimizing interference from other signals. An LFF of 10 Hz helps remove very low-frequency noise and movement artifacts that are not related to muscle activity but could potentially distort the EMG readings. Since muscle activity generates signals within a range that starts from around 10 Hz, this setting allows for adequate detection of the EMG signal without losing important data. Setting the HFF at 100 Hz ensures that higher frequency noise, which may arise from electrical interference or rapid muscle fiber firing, is filtered out, allowing for a clearer representation of the actual muscle activity. Frequencies above 100 Hz are typically less representative of the muscle contraction signals of interest, so this high cutoff effectively filters out irrelevant data without impacting the critical information needed to analyze muscle function. Using these specific filter settings supports the accurate interpretation of EMG signals essential for diagnosing muscular and neurological conditions, which is vital in the practice of sleep technology and related fields.*

SAMPLE

# Next Steps

Congratulations on reaching the final section of this guide. You've taken a meaningful step toward passing your certification exam and advancing your career.

As you continue preparing, remember that consistent practice, review, and self-reflection are key to success. Make time to revisit difficult topics, simulate exam conditions, and track your progress along the way.

If you need help, have suggestions, or want to share feedback, we'd love to hear from you. Reach out to our team at hello@examzify.com.

Or visit your dedicated course page for more study tools and resources:

https://registeredsleeptechnologist.examzify.com

We wish you the very best on your exam journey. You've got this!

SAMPLE